TinyScale

Create happy beginnings

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Baby weight

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User manual

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Security:

Before using the product, read this user manual and follow the recommendations contained therein when using the device.

Device description [page 2 - A]:

1. On/Off/Tare, 2. Unit conversion, 3. LCD screen, 4. Weight retention.

Features [page 2 - A]:

1. ON/OFF button – press to turn the scale on or off and then press again to activate the TARE function.
2. TARE button – press to zero the scale before placing clothes and blankets on the tray. This will allow you to subtract their weight and get your baby's exact weight.
3. When "0" appears on the display, you can place your baby on the scale.
4. HOLD button – once the weight value has stabilized, press the button to save the weight on the screen. To cancel saved data, press the button again.
5. UNIT button – allows you to change mass units with one press.
6. Auto-off function (energy-saving mode) - the device turns off automatically after about 2 minutes of inactivity.
7. Low Battery Warning - When "LO" appears on the screen, replace the battery.

How to use the scale:

1. Place the scale on a hard, flat surface.
2. Press the on/off button to turn it on, then press "TARE" so that the display shows "0".
3. When the display shows "0", place your baby with his full weight on the scale. Press "UNIT" button to select unit change (g - kg - oz - lb:oz).

Warnings:

1. Before use, place the scale on a flat surface.
2. Do not place or use the product in humid, very hot or very cold places.
3. Protect the product from damage due to falling or strong shock as it is a high precision electronic product.
4. To clean the product, use only a soft cloth, without water or other chemical solutions.
5. This product is intended for measurements in a home environment and cannot be used for professional measurements.
6. Remove the battery if you will not use the product for a long time.
7. Remove items from the scale to avoid damage if ERR appears on the display (maximum weight is 30 kg).
8. Do not leave your child alone on the scale without supervision.
9. The product is not a toy, keep it out of the reach of children.

Power supply: 2 AAA batteries

Measuring range: 10g - 30000g
Precision: 2g
Measurement deviation: $\pm 0,3\%$
Product weight: 1.3kg
Length measurement: 0 cm - 60 cm

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[EN] Rules for handling electrical and electronic equipment: Waste electrical and electronic equipment as well as batteries and accumulators cannot be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage may result in health damage or environmental pollution. Proper household handling of waste electrical and electronic equipment and batteries and accumulators positively contributes to the reuse and recovery, including recycling, of waste equipment. Waste electrical and electronic equipment as well as batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and performed the same functions as the equipment sold. When delivering equipment intended for households to the buyer, the distributor is obliged to collect used equipment from households free of charge at the place of delivery of the equipment, provided that the used equipment is of the same type and performed the same functions as the delivered equipment. A distributor running a retail unit with a sales area within the meaning of Art. 2 point 19 of the Act of March 27, 2003 on spatial planning and development of at least 400 m² devoted to the sale of equipment intended for households, is obliged to accept waste equipment from households, free of charge, in this unit or in its immediate vicinity, none of the external dimensions exceeds 25 cm, without the need to purchase new household equipment.

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[EN] Summary declaration of conformity: The manufacturer declares, under its sole responsibility, that the device complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances. The full content of the declaration of conformity is available at www.momini.eu/deklaracje

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EN | More information about the product is available on the manufacturer's website at:

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