BEAUTIFLY

SkinGlow

PRO

series

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Rated power: 18 W
Actual operating power: ≤ 4.5 W
Power supply: DC 5 V / 1 A
Device dimensions: 27.5 × 20 × 3.5 cm
Weight: approx. 1.37 kg

Safety:

Before using the product, read this instruction manual carefully and follow the recommendations during operation.

Package contents:

1. LED phototherapy panel / 2. Protective goggles (1 pc.) / 3. USB Type-C power cable (1 pc.) / 4. Instruction manual

Device structure [page 2 - f1]:

1. LED diode matrix with different wavelengths / 2. Digital display / 3. Time + increase operating time / 4. Time – decrease operating time / 5. M button mode / light color change from M1 to M12 / 6. ON/OFF button device power on and off

Device operation

1. Place the panel on a stable surface / 2. Connect the device to power / 3. Press and hold the ON/OFF button for about 1.5 seconds to turn on the device / 4. Press the M button and select one of the modes from M1 to M12 to change the color / light mode (Mode M12 activates automatic color cycling) / 5. Set the working time: press Time – to shorten the treatment time and Time + to extend it

NOTE: After the countdown ends, the device will switch off automatically. To turn off the device manually, hold the ON/OFF button for approximately 1.5–2 seconds.

Method of use:

1. Clean the skin before treatment / 2. You may apply a suitable skincare product (without strong photosensitizing ingredients) / 3. Put on the included protective goggles / 4. Adjust the appropriate panel tilt angle [page 2 - f2] / 5. Maintain a safe distance from the light source / 6. Turn on the device and begin the session.

Charging, maintenance and care:

If the device does not respond or the battery level is low, connect the cable to the USB Type-C port and a 5 V / 1 A power source, and disconnect after full charging. Avoid direct contact with water. Clean with a dry or slightly damp cloth. Store in a dry place away from direct sunlight.

Precautions:

Keep out of reach of children. The device is not intended for persons under 18 years of age, pregnant women, women planning pregnancy or breastfeeding, individuals sensitive to light, or persons with genetic or skin diseases (e.g. lupus, eczema, albinism). Do not use while charging. Do not drop, hit, or apply excessive pressure. Do not disassemble or modify the device. Do not immerse in water or clean with corrosive agents. Store in a cool, dry place. Do not use on open wounds, rashes, swollen, irritated, infected or inflamed skin, or on areas with increased risk of skin cancer. When using IPL or laser devices, wear eye protection. For 48 hours before and after treatment, avoid tanning beds and intense sun exposure. Redness or irritation may occur. If this happens, discontinue use and consult a doctor. More frequent use than recommended increases the risk of irritation. Before starting therapy, medical consultation and, if necessary, diagnostic tests are recommended. Do not perform the procedure without documentation from such consultation.

More information about the company and products [page 2 - f3]:

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[EN] Rules for handling electrical and electronic equipment: Used electrical and electronic equipment, batteries, and accumulators should not be disposed of with other waste. This is indicated by the crossed-out waste bin symbol on the product. Electrical and electronic equipment, batteries, and accumulators may contain accessories, components, and parts that are harmful to the environment and human health. Improper handling or their consequences may result in damage to the environment or the environment. Proper household handling of used electrical and electronic equipment, batteries, and accumulators is essential for the proper disposal of devices and equipment, including used electrical and electronic equipment. Used electrical and electronic equipment, batteries, and accumulators can be disposed of separately. The distributor is obligated to collect household equipment free of charge at the point of sale, provided the used equipment is of the same type and performs the same functions as commercial equipment. A distributor delivering household appliances to a purchaser is obligated to collect waste household appliances free of charge at the point of delivery, provided the waste appliances are of the same type and have served the same function as the appliances. A distributor operating a retail business with an online sales area of 400 m² of household appliances designated for sale pursuant to Article 2, Section 19 of the Act of 27 March 2003 on Spatial Planning and Development is entitled to free membership in the same entity or its immediate vicinity for discarded household appliances with external dimensions not exceeding 25 cm, without the need to purchase new household appliances. The product was registered on the market after August 13, 2005.

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[EN] Summary declaration of conformity: The manufacturer declares on its sole responsibility that the product listed below complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility, standards used: EN IEC 55014-1:2021, EN IEC 55014-2:2021, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances and Directive (EU) 2017/2102 of the European Parliament and of the Council of 15 November 2017 amending Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment, applied standards: IEC62321-3-1:2013, IEC62321-4:2013+A1:2017, IEC62321-5:2013, IEC62321-6:2015, IEC62321-7-1:2015, IEC62321-7-2:2017, IEC62321-8:2017. The full text of the declaration of conformity is available on the website www.beautifly.eu/deklaracje

Importer: Beautifly Sp. z o.o. Plac Bankowy 2, 00-095 Warszawa, Polska, kontakt@beautifly.eu
Model: SkinGlow Pro
SN: 02265903003700913
Marka: Beautifly
e-mail: serwis@beautifly.eu
Wyprodukowano w Chinach